




ROTOR SPECIFICATIONS

Type	Fixed Angle
Material	Aluminum
Capacity	6 x 500 mL
Maximum Speed	12,000 rpm
Tube Dimensions Ø x L	70 x 190 mm
Angle	23°
K Factor	2364
Net Weight	14.4 kg/31.75 lbs
Warranty	10 years
Certifications¹	

	RCF	RADIUS
Maximum	24,336 xg	151.3 mm
Minimum	6,337 xg	39.4 mm

ROTOR AND ACCESSORIES ORDERING INFORMATION


		SORVALL RC-6 Plus Series	THERMO SCIENTIFIC SORVALL Evolution RC Series	SORVALL RC-5 Plus Series
SLA-3000	Max. Speed	11,000 rpm	12,000 rpm	11,000 rpm
	Max. RCF	20,450 xg	24,340 xg	20,450 xg

SLA-3000 Rotor Package includes:

Description	Cat. No.
SLA-3000 Rotor	07149
Rotor Cover Assembly	24512
Vacuum Grease	76003500
Lubricant	70006692
Rotor Extraction Tool	01047
Instruction Manual	11463
Rotor Cover O-ring (qty. 6)	92232
Retaining Ring Pliers	65614
Locking Knob O-ring (qty. 6)	90442
Center Hold O-ring (qty. 6)	92231

LABWARE ORDERING INFORMATION

Tube Vol. (mL)	Vol. Fill (mL)	Description	Cat. No.	Qty. /Pkg.	Tubes /Rotor	Max. Speed (rpm)	Max. RCF (xg)	Dims ØxL (mm)	REQUIRED TUBE CLOSURES, ADAPTERS, TOOLS & ACCESSORIES			
									Type	Cat. No.	Qty. /Pkg.	Description
500.0	450.0	Nalgene® PC Oak Ridge Bottle	3140-0500	4	6	9,000	13,689	70 x 160	Closure	Included	4	PP Sealing
500.0	450.0	Nalgene PP Oak Ridge Bottle	3141-0500	4	6	9,000	13,689	70 x 160	Closure	Included	4	PP Sealing

¹  Biocontainment certification by the Public Health Laboratory Service, Centre for Applied Microbiological Research, Porton Down, UK, shown to be suitable for use with materials up to ACDP Category 3 as categorized by the Advisory Committee on Dangerous Pathogens.

LABWARE ORDERING INFORMATION

Tube Vol. (mL)	Vol. Fill (mL)	Description	Cat. No.	Qty. /Pkg.	Tubes /Rotor	Max. Speed (rpm)	Max. RCF (xg)	Dims ØxL (mm)	REQUIRED TUBE CLOSURES, ADAPTERS, TOOLS & ACCESSORIES			
									Type	Cat. No.	Qty. /Pkg.	Description
500.0	480.0	PP Dry-Spin Bottle	06829	6	6	9,000	13,689	70 x 190	Closure	Included	6	PP Sealing
									Closure ²	03545	24	Replacement O-rings
									Closure ²	06841	12	Replacement PP Sealing Lids
250.0	270.0	PC Dry-Spin Bottle	03098	12	6	12,000	24,336	61 x 146	Closure	Included	12	PP Sealing
									Closure ²	03545	24	Replacement O-rings
									Closure ²	06841	12	Replacement PP Sealing Lids
									Adapter	41002	1	1 place/adapter
250.0	280.0	PP Dry-Spin Bottle	03069	12	6	10,000	16,900	61 x 146	Closure	Included	12	PP Sealing
									Closure ²	03545	24	Replacement O-rings
									Closure ²	06841	12	Replacement PP Sealing Lids
									Adapter	41002	1	1 place/adapter
250.0	250.0	Nalgene PC Oak Ridge Bottle	3140-0250	4	6	12,000	24,336	61 x 135	Closure	Included	4	PP Sealing
									Adapter	41002	1	1 places/adapter
250.0	250.0	Nalgene PP Oak Ridge Bottle	3141-0250	4	6	12,000	24,336	61 x 133	Closure	Included	4	PP Sealing
									Adapter	41002	1	1 place/adapter
100.0	100.0	GL Pyrex [®] Tube	03230	1	6	3,000	1,521	44 x 137	Closure	Included	1	PP Non-Sealing
									Adapter	00371 ³	1	1 place adapter
									Adapter	41002 ³	1	1 place adapter
50.0	50.0	Nunc [™] PP Conical Tissue Culture Tube	334959	25	6	5,450	5,020	30 x 55	Closure	Included	25	HDPE Sealing
									Adapter	03072	1	1 place/adapter
									Adapter	410022	1	1 place/adapter

² Replacement parts. Not required for use.

³ Both adapters must be used together.

© 2010 Thermo Fisher Scientific Inc. All rights reserved. Pyrex is a registered trademark of Corning Inc. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

For more information, please visit
www.thermoscientific.com/rotors.

Thermo
 SCIENTIFIC